

THE COLEMAN™ MICROINJECTION SYSTEM



body aesthetics

THE COLEMAN™ MICROINJECTION SYSTEM

The COLEMAN™ Microinjection System is an innovative concept in instrumentation that assists with microinfiltration procedures including donor site preparation, removal and preparation of aspirant, and precise infiltration at the recipient site.

Designed by Dr. Sydney Coleman and perfected in clinical use, this instrumentation offers everything you need for small volume fat transfer procedures.

CANNULA

COLEMAN™ Infiltration Cannula



Three styles and two sizes have been developed to perform soft tissue infiltration.

COLEMAN™ Aspiration Cannula



Configured to allow passage of aspirant optimally sized for each infiltration procedure.

COLEMAN™ “V” Dissector Cannula



An atraumatic “V” channel at the distal tip helps disrupt soft tissue for controlled release of adhesions and ligamentous connections.


COLEMAN™ Contour Series



Designed to access difficult deep or superficial graft sites, these concave and convex shaft styles are available with Style I, II and III tips.

COLEMAN™ Style Mini Infiltration Cannula

The COLEMAN™ Style Mini Infiltration Cannula is ideal for areas requiring precise control. Available in two sizes, I (19 ga) and II (20 ga).



TIP STYLES

Style I

COLEMAN™ Style I features a blunt closed end distally, with a single opening.

Style II

COLEMAN™ Style II features a blunt end, but is more open distally and laterally than Style I. The Style II is the cannula most frequently used by Dr. Coleman.

Style III

COLEMAN™ Style III is open and “flattened” distally, and is used by Dr. Coleman when working superficially.

Aspiration

COLEMAN™ Aspiration Cannula allows passage of aspirant optimally sized for each infiltration procedure.

“V” Dissector

Coleman™ “V” Dissector Cannula has an atraumatic “V” channel at the distal tip.

ACCESSORIES


COLEMAN™ Syringe Caps

Positive connections assure a leak proof fit to luer lock syringes during the refinement process. Caps are designed to withstand the fluid pressures generated during centrifugation.



COLEMAN™ Syringe Rack

Designed to hold 1cc, 3cc and 10cc syringes, this autoclavable rack is used by Dr. Coleman to facilitate preparation and counting during his procedure.



MediLite® Centrifuge

The MediLite® Centrifuge is capable of spinning six 10cc syringes simultaneously at 3100 rpm. Rotor and tubes are autoclavable.



*The third-party trademarks used herein are trademarks of their respective owners.

COLEMAN™ PRODUCTS ORDERING INFORMATION

Catalog #	Description	Length	Quantity
COLEMAN™ INFILTRATION CANNULA			
COL-I7	COLEMAN™ Infiltration Cannula, Style I	7cm	EA
COL-II7	COLEMAN™ Infiltration Cannula, Style II	7cm	EA
COL-III7	COLEMAN™ Infiltration Cannula, Style III	7cm	EA
COL-II7C	COLEMAN™ Curved Infiltration Cannula, Style II	7cm	EA
COL-I9	COLEMAN™ Infiltration Cannula, Style I	9cm	EA
COL-II9	COLEMAN™ Infiltration Cannula, Style II	9cm	EA
COL-III9	COLEMAN™ Infiltration Cannula, Style III	9cm	EA
COL-I14B	COLEMAN™ Infiltration Cannula, Style I for Body - 14g	15cm	EA
COL-I16B	COLEMAN™ Infiltration Cannula, Style I for Body - 16g	15cm	EA
COL-I16BC	COLEMAN™ Curved Infiltration Cannula, Style I for Body - 16g	15cm	EA
COL-STYLET	COLEMAN™ Stylet		EA
COLEMAN™ ASPIRATION CANNULA			
COL-ASP I	Aspiration Cannula	15cm	EA
COL-ASP IC	Curved Aspiration Cannula	15cm	EA
COL-ASP B	Aspiration Cannula, Body	23cm	EA
COL-ASP BC	Curved Aspiration Cannula, Body	23cm	EA
COLEMAN™ “V” DISSECTOR CANNULA			
COL-V	COLEMAN™ “V” Dissector		10/box
COL-V149	COLEMAN™ “V” Dissector - 14g	9cm	15/box
COL-V12B	COLEMAN™ “V” Dissector - 12g	15cm	10/box
COL-V14B	COLEMAN™ “V” Dissector - 14g	15cm	10/box
COLEMAN™ CONTOUR SERIES CANNULA			
COL-I9CC	COLEMAN™ Concave Infiltration Cannula, Style I	9cm	EA
COL-II9CC	COLEMAN™ Concave Infiltration Cannula, Style II	9cm	EA
COL-III9CC	COLEMAN™ Concave Infiltration Cannula, Style III	9cm	EA
COL-I9CV	COLEMAN™ Convex Infiltration Cannula, Style I	9cm	EA
COL-II9CV	COLEMAN™ Convex Infiltration Cannula, Style II	9cm	EA
COL-III9CV	COLEMAN™ Convex Infiltration Cannula, Style III	9cm	EA
COL-II9S	COLEMAN™ “S” Infiltration Cannula, Style II	9cm	EA
COLEMAN™ STYLE MINI INFILTRATION CANNULA			
COL-MINI I-S	COLEMAN™ Mini I Stylet		EA
COL-MINI I-K	COLEMAN™ Mini I Kit		EA
COL-MINI II-S	COLEMAN™ Mini II Stylet		EA
COL-MINI II-K	COLEMAN™ Mini II Kit		EA
COLEMAN™ ACCESSORIES			
COL-1223C	COLEMAN™ Curved Irrigation Cannula - 12g X 23cm		EA
COL-LLCAP	COLEMAN™ Luer Lock Caps, Sterile		6/Pkg, 25 Pkg/Bx
COL-RACK	Autoclavable Syringe Rack		EA
B9628	1cc BD Luer Lock Syringe, Sterile		100/Bx
B9604	10cc BD Luer Lock Syringe, Sterile		100/Bx
CENTRIFUGE			
2580	MediLite® Centrifuge, 110 V		EA
2581	MediLite® Centrifuge, 220 V		EA
MS-TUBE	Extra Tube for MediLite® Centrifuge		EA
MS-ROTOR	Extra Rotor for MediLite® Centrifuge		EA

*The third-party trademarks used herein are trademarks of their respective owners.



201 Mentor Drive
Santa Barbara, CA 93111 USA
Tel: 805-879-6000
www.mentorworldwide.com

Customer Service
Phone: 800-235-5731
Fax: 805-967-7108