

Consistency in Excellence.

Contour Profile® Gel Family – The 300 Series.



CPG™ 331



CPG™ 332



CPG™ 333



CPG™ 321



CPG™ 322



CPG™ 323



CPG™ 311



CPG™ 312



CPG™ 313



CPG™ 300 series. Simplicity itself.

CPG™ NAMING LEGEND

Cohesive Level	Implant Height	Implant Projection
III 3	Tall 3	High 3
	Medium 2	Moderate <i>Plus</i> 2
	Low 1	Moderate 1
KEY: Digit 1= Cohesive Level Digit 2= Height of Implant Digit 3= Projection of Implant		

Introducing a simple way to identify exactly what cohesive level, height and projection each implant provides. Three numbers are all you need to know.

The first number identifies the implant type. Contour Profile® is the 300 series, and an easy way to remember that all CPG™ implants are Cohesive III™. The second number relates to the height of the implant, with 1 being the lowest height and 3 being the tallest. The third number identifies implant projection, with 1 representing Moderate, 2 representing Moderate *Plus* and 3 representing High Projection. Ordering is as simple as 1 – 2 – 3.

A smarter approach to implant selection.

Implant selection should be as simple as it is accurate. With that goal in mind, Mentor developed The Contour Profile® Gel (CPG™) 300 Series. Three heights, three projections and multiple ways to match your patient's unique physiology.

Designed to fit.

Contour Profile® Gel implants combine a streamlined inventory with highly flexible sizing. Essentially, physicians require less stock to match more body types.

- Three Projections: Moderate, Moderate *Plus* and High
- Three Heights: Low, Medium and Tall
- Unique shape: A projection point designed to maximize visual outcome
- Cohesive III™ MemoryGel®, Mentor's most cohesive, form stabilizing gel
- Gradual upper pole slope for implant stability
- Orientation marks to facilitate surgical placement
- Siltex™ texturing



A cross-section of a Contour Profile® implant, Cohesive III™

Mentor's proprietary cohesive silicone gel holds shape uniformly and resembles the feel of actual breast tissue. Beautifully pliant. Firm without rigidity. MemoryGel® offers patients a more natural look and feel. By varying gel components, Mentor produces a wide selection of consistencies from soft (Cohesive I™) to firm (Cohesive III™).

All CPG™ products feature Cohesive III™ MemoryGel®, providing the optimal consistency for aesthetic results.

Important Safety Information:

MemoryGel® breast implants are indicated for breast augmentation in women at least 18 years old or for breast reconstruction. Breast implant surgery should not be performed in women with active infection anywhere in their body with existing cancer or pre-cancer of their breast who have not received adequate treatment for those conditions or are pregnant or nursing.

Breast implants are not lifetime devices and breast implantation is not necessarily a one-time surgery. The most common complications with the MemoryGel® breast implants include reoperation, capsular contracture, asymmetry, and breast pain. A lower risk of complication is rupture. The health consequences of a ruptured silicone gel-filled breast implant have not been fully established. Screenings such as mammography, MRI, or ultrasound are recommended after initial implant surgery to assist in detecting implant rupture.


Your patient needs to be informed and understand the risks and benefits of breast implants, and provided with an opportunity to consult with you prior to deciding on surgery.

For detailed indications, contraindications, warning and precautions associated with the use of MemoryGel® breast implants please refer to the Product Insert Data Sheet provided with each product, or review the Important Safety Information provided at www.mentorcorp.com.

Contour Profile® Gel Family – The 300 Series.

Featuring Cohesive III™ MemoryGel®.

CPG™ NAMING LEGEND		
Cohesive Level	Implant Height	Implant Projection
III 3	Tall 3	High 3
	Medium 2	Moderate <i>Plus</i> 2
	Low 1	Moderate 1
KEY: Digit 1 = Cohesive Level Digit 2 = Height of Implant Digit 3 = Projection of Implant		



CPG™ 331

Cohesive III™, Tall Height,
Moderate Projection

The Tall Height, Moderate projection CPG™ implant is designed for tall patients with a narrow torso seeking for moderate results in their augmentation procedure. This implant also shows excellent results when correcting tuberous breasts in the above mentioned patient type.

As a reconstruction device it can be used as a follow up implant after the Tall Height Contour Profile® Breast Tissue Expander.

CPG™ 332

Cohesive III™, Tall Height,
Moderate *Plus* Projection

The Tall Height, Moderate Plus projection CPG™ implant has the same footprint as the CPG™ 331 but with more inferior pole projection. It is an ideal device for patients with taller, narrow torsos desiring more projection.

This implant is also selected for reconstruction patients after using Tall Height Contour Profile® Tissue Expanders.

CPG™ 333

Cohesive III™, Tall Height,
High Projection

The Tall Height, High projection CPG™ implant provides tall height with the most projection. It is designed for augmentation procedures in tall, slender patients with a narrow torso when a high projection outcome is desired.

In a two stage reconstruction procedure the implant can be used after our Tall Height Contour Profile® Breast Tissue Expanders. It is the ideal choice as well when having to match a large contralateral breast in single stage reconstruction cases.

CPG™ 321

Cohesive III™, Medium Height,
Moderate Projection

The Medium Height, Moderate projection CPG™ implant is an excellent choice for most augmentation patients whether the augmentation procedure is associated with mastopexy or not.

It also provides an option in reconstruction patients when a high projection is not desired and can be used as a follow up implant after our Medium Height Contour Profile® Breast Tissue Expander.

CPG™ 322

Cohesive III™, Medium Height,
Moderate *Plus* Projection

The Medium Height, Moderate Plus projection CPG™ implant is suitable for a great number of augmentation patients. It offers Moderate Plus projection especially in the inferior pole.

As a reconstruction device it can be used after our Medium Height Contour Profile® Breast Tissue Expanders or associated with a Lat. Dorsi flap.

CPG™ 323

Cohesive III™, Medium Height,
High Projection

The Medium Height, High projection CPG™ implant offers multiple options in augmentation procedures especially when correcting breast ptosis and no mastopexy is associated.

As well, it is an ideal device for reconstruction procedures where a high inferior pole projection is important. In a two stage reconstruction procedure the implant can be used after our Medium Height Contour Profile® Breast Tissue Expanders.

CPG™ 311

Cohesive III™, Low Height,
Moderate Projection

The Low Height, Moderate projection CPG™ implant is designed for augmentation or reconstruction procedures in short patients with a wide torso seeking regular projection.

As a reconstruction device it can be used after both our Low Height and Medium Height Contour Profile® Breast Tissue Expanders.

CPG™ 312

Cohesive III™, Low Height,
Moderate *Plus* Projection

The Low Height, Moderate Plus projection CPG™ implant with greater inferior pole projection offers multiple options in breast surgery. It is a highly suitable implant for patients with a shorter, wider torso seeking breast augmentation with Moderate Plus projection especially in the lower pole.

As a reconstruction device it can be used after both our Low Height and Medium Height Contour Profile® Breast Tissue Expanders.

CPG™ 313

Cohesive III™, Low Height,
High Projection

The Low Height, High projection CPG™ implant is an ideal augmentation device for short patients with a wide torso seeking high inferior pole projection.

In a two stage reconstruction procedure the implant can be used after both our Low Height and Medium Height Contour Profile® Breast Tissue Expanders.

Contour Profile® Gel Family – The 300 Series



CPG™ 331

Cohesive III™, Tall Height, Moderate Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
125 cc	9,0 cm	9,2 cm	3,2 cm	334-0903	RSZ-0903
150 cc	9,5 cm	9,7 cm	3,3 cm	334-0953	RSZ-0953
175 cc	10,0 cm	10,2 cm	3,4 cm	334-1003	RSZ-1003
200 cc	10,5 cm	10,7 cm	3,5 cm	334-1053	RSZ-1053
230 cc	11,0 cm	11,3 cm	3,6 cm	334-1103	RSZ-1103
265 cc	11,5 cm	11,8 cm	3,7 cm	334-1153	RSZ-1153
300 cc	12,0 cm	12,3 cm	3,9 cm	334-1203	RSZ-1203
340 cc	12,5 cm	12,8 cm	4,0 cm	334-1253	RSZ-1253
380 cc	13,0 cm	13,3 cm	4,1 cm	334-1303	RSZ-1303
425 cc	13,5 cm	13,8 cm	4,3 cm	334-1353	RSZ-1353
475 cc	14,0 cm	14,3 cm	4,5 cm	334-1403	RSZ-1403
530 cc	14,5 cm	14,8 cm	4,6 cm	334-1453	RSZ-1453
585 cc	15,0 cm	15,3 cm	4,8 cm	334-1503	RSZ-1503
645 cc	15,5 cm	15,9 cm	5,0 cm	334-1553	RSZ-1553

CPG™ 321

Cohesive III™, Medium Height, Moderate Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
120 cc	9,0 cm	8,5 cm	3,3 cm	354-0908	RSZ-0908
135 cc	9,5 cm	8,9 cm	3,5 cm	354-0958	RSZ-0958
155 cc	10,0 cm	9,4 cm	3,7 cm	354-1008	RSZ-1008
180 cc	10,5 cm	9,9 cm	3,8 cm	354-1058	RSZ-1058
215 cc	11,0 cm	10,3 cm	3,9 cm	354-1108	RSZ-1108
245 cc	11,5 cm	10,8 cm	4,0 cm	354-1158	RSZ-1158
280 cc	12,0 cm	11,3 cm	4,2 cm	354-1208	RSZ-1208
315 cc	12,5 cm	11,8 cm	4,4 cm	354-1258	RSZ-1258
355 cc	13,0 cm	12,2 cm	4,6 cm	354-1308	RSZ-1308
395 cc	13,5 cm	12,7 cm	4,7 cm	354-1358	RSZ-1358
440 cc	14,0 cm	13,2 cm	4,9 cm	354-1408	RSZ-1408
480 cc	14,5 cm	13,6 cm	5,0 cm	354-1458	RSZ-1458
530 cc	15,0 cm	14,1 cm	5,2 cm	354-1508	RSZ-1508
640 cc	16,0 cm	15,0 cm	5,6 cm	354-1608	RSZ-1608
775 cc	17,0 cm	16,0 cm	5,9 cm	354-1708	RSZ-1708

CPG™ 332

Cohesive III™, Tall Height, Moderate *Plus* Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
145 cc	9,0 cm	9,4 cm	3,8 cm	334-0909	RSZ-0909
175 cc	9,5 cm	9,9 cm	4,0 cm	334-0959	RSZ-0959
205 cc	10,0 cm	10,4 cm	4,2 cm	334-1009	RSZ-1009
235 cc	10,5 cm	10,9 cm	4,4 cm	334-1059	RSZ-1059
270 cc	11,0 cm	11,5 cm	4,7 cm	334-1109	RSZ-1109
305 cc	11,5 cm	12,0 cm	4,9 cm	334-1159	RSZ-1159
350 cc	12,0 cm	12,5 cm	5,1 cm	334-1209	RSZ-1209
395 cc	12,5 cm	13,0 cm	5,3 cm	334-1259	RSZ-1259
445 cc	13,0 cm	13,5 cm	5,5 cm	334-1309	RSZ-1309
495 cc	13,5 cm	14,1 cm	5,7 cm	334-1359	RSZ-1359
555 cc	14,0 cm	14,6 cm	5,9 cm	334-1409	RSZ-1409
615 cc	14,5 cm	15,1 cm	6,1 cm	334-1459	RSZ-1459
680 cc	15,0 cm	15,6 cm	6,3 cm	334-1509	RSZ-1509

CPG™ 322

Cohesive III™, Medium Height, Moderate *Plus* Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
140 cc	9,0 cm	8,5 cm	3,8 cm	334-0905	RSZ-0905
165 cc	9,5 cm	8,9 cm	4,0 cm	334-0955	RSZ-0955
195 cc	10,0 cm	9,4 cm	4,2 cm	334-1005	RSZ-1005
225 cc	10,5 cm	9,9 cm	4,4 cm	334-1055	RSZ-1055
255 cc	11,0 cm	10,3 cm	4,7 cm	334-1105	RSZ-1105
295 cc	11,5 cm	10,8 cm	4,9 cm	334-1155	RSZ-1155
330 cc	12,0 cm	11,3 cm	5,1 cm	334-1205	RSZ-1205
375 cc	12,5 cm	11,8 cm	5,3 cm	334-1255	RSZ-1255
420 cc	13,0 cm	12,2 cm	5,5 cm	334-1305	RSZ-1305
475 cc	13,5 cm	12,7 cm	5,7 cm	334-1355	RSZ-1355
525 cc	14,0 cm	13,2 cm	5,9 cm	334-1405	RSZ-1405
585 cc	14,5 cm	13,6 cm	6,1 cm	334-1455	RSZ-1455
650 cc	15,0 cm	14,1 cm	6,3 cm	334-1505	RSZ-1505

CPG™ 333

Cohesive III™, Tall Height, High Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
180 cc	9,0 cm	9,4 cm	4,6 cm	334-0904	RSZ-4090
215 cc	9,5 cm	9,9 cm	4,8 cm	334-0954	RSZ-4095
250 cc	10,0 cm	10,4 cm	5,1 cm	334-1004	RSZ-4100
290 cc	10,5 cm	10,9 cm	5,3 cm	334-1054	RSZ-4105
330 cc	11,0 cm	11,5 cm	5,6 cm	334-1104	RSZ-4110
380 cc	11,5 cm	12,0 cm	5,8 cm	334-1154	RSZ-4115
430 cc	12,0 cm	12,5 cm	6,0 cm	334-1204	RSZ-4120
485 cc	12,5 cm	13,0 cm	6,2 cm	334-1254	RSZ-4125
545 cc	13,0 cm	13,5 cm	6,5 cm	334-1304	RSZ-4130
610 cc	13,5 cm	14,1 cm	6,7 cm	334-1354	RSZ-4135
680 cc	14,0 cm	14,6 cm	6,9 cm	334-1404	RSZ-4140
755 cc	14,5 cm	15,1 cm	7,1 cm	334-1454	RSZ-4145

CPG™ 323

Cohesive III™, Medium Height, High Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
165 cc	9,0 cm	8,5 cm	4,6 cm	334-0902	RSZ-0902
195 cc	9,5 cm	8,9 cm	4,8 cm	334-0952	RSZ-0952
225 cc	10,0 cm	9,4 cm	5,1 cm	334-1002	RSZ-1002
260 cc	10,5 cm	9,9 cm	5,3 cm	334-1052	RSZ-1052
300 cc	11,0 cm	10,3 cm	5,6 cm	334-1102	RSZ-1102
345 cc	11,5 cm	10,8 cm	5,8 cm	334-1152	RSZ-1152
390 cc	12,0 cm	11,3 cm	6,0 cm	334-1202	RSZ-1202
440 cc	12,5 cm	11,8 cm	6,2 cm	334-1252	RSZ-1252
495 cc	13,0 cm	12,2 cm	6,5 cm	334-1302	RSZ-1302
555 cc	13,5 cm	12,7 cm	6,7 cm	334-1352	RSZ-1352
620 cc	14,0 cm	13,2 cm	6,9 cm	334-1402	RSZ-1402
685 cc	14,5 cm	13,6 cm	7,1 cm	334-1452	RSZ-1452



CPG™ 311

Cohesive III™, Low Height, Moderate Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
100 cc	9,0 cm	8,0 cm	3,0 cm	334-0901	RSZ-0090
120 cc	9,5 cm	8,5 cm	3,2 cm	334-0951	RSZ-0095
140 cc	10,0 cm	8,9 cm	3,4 cm	334-1001	RSZ-0100
160 cc	10,5 cm	9,4 cm	3,6 cm	334-1051	RSZ-0105
180 cc	11,0 cm	9,8 cm	3,7 cm	334-1101	RSZ-0110
210 cc	11,5 cm	10,3 cm	3,9 cm	334-1151	RSZ-0115
235 cc	12,0 cm	10,7 cm	4,1 cm	334-1201	RSZ-0120
270 cc	12,5 cm	11,2 cm	4,2 cm	334-1251	RSZ-0125
300 cc	13,0 cm	11,6 cm	4,4 cm	334-1301	RSZ-0130
335 cc	13,5 cm	12,0 cm	4,6 cm	334-1351	RSZ-0135
375 cc	14,0 cm	12,5 cm	4,7 cm	334-1401	RSZ-0140
415 cc	14,5 cm	12,9 cm	4,9 cm	334-1451	RSZ-0145
460 cc	15,0 cm	13,4 cm	5,1 cm	334-1501	RSZ-0150
510 cc	15,5 cm	13,8 cm	5,2 cm	334-1551	RSZ-0155
560 cc	16,0 cm	14,3 cm	5,4 cm	334-1601	RSZ-0160
615 cc	16,5 cm	14,7 cm	5,6 cm	334-1651	RSZ-0165

CPG™ 312

Cohesive III™, Low Height, Moderate *Plus* Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
125 cc	9,0 cm	8,0 cm	3,8 cm	334-0907	RSZ-0907
145 cc	9,5 cm	8,4 cm	4,0 cm	334-0957	RSZ-0957
170 cc	10,0 cm	8,8 cm	4,2 cm	334-1007	RSZ-1007
195 cc	10,5 cm	9,3 cm	4,4 cm	334-1057	RSZ-1057
225 cc	11,0 cm	9,7 cm	4,7 cm	334-1107	RSZ-1107
255 cc	11,5 cm	10,2 cm	4,9 cm	334-1157	RSZ-1157
290 cc	12,0 cm	10,6 cm	5,1 cm	334-1207	RSZ-1207
330 cc	12,5 cm	11,1 cm	5,3 cm	334-1257	RSZ-1257
370 cc	13,0 cm	11,5 cm	5,5 cm	334-1307	RSZ-1307
415 cc	13,5 cm	12,0 cm	5,7 cm	334-1357	RSZ-1357
465 cc	14,0 cm	12,4 cm	5,9 cm	334-1407	RSZ-1407
515 cc	14,5 cm	12,8 cm	6,1 cm	334-1457	RSZ-1457
570 cc	15,0 cm	13,3 cm	6,3 cm	334-1507	RSZ-1507
690 cc	16,0 cm	14,2 cm	6,8 cm	334-1607	RSZ-1607

CPG™ 313

Cohesive III™, Low Height, High Projection

Vol.	Width	Height	Proj.	Catalog #	Resterilizable Gel Sizer*
130 cc	9,0 cm	8,1 cm	4,4 cm	334-0906	RSZ-0906
155 cc	9,5 cm	8,6 cm	4,5 cm	334-0956	RSZ-0956
180 cc	10,0 cm	9,0 cm	4,6 cm	334-1006	RSZ-1006
210 cc	10,5 cm	9,5 cm	4,8 cm	334-1056	RSZ-1056
240 cc	11,0 cm	9,9 cm	4,9 cm	334-1106	RSZ-1106
270 cc	11,5 cm	10,3 cm	5,0 cm	334-1156	RSZ-1156
310 cc	12,0 cm	10,8 cm	5,2 cm	334-1206	RSZ-1206
350 cc	12,5 cm	11,2 cm	5,4 cm	334-1256	RSZ-1256
395 cc	13,0 cm	11,7 cm	5,6 cm	334-1306	RSZ-1306
440 cc	13,5 cm	12,2 cm	5,8 cm	334-1356	RSZ-1356
490 cc	14,0 cm	12,6 cm	6,1 cm	334-1406	RSZ-1406
545 cc	14,5 cm	13,1 cm	6,3 cm	334-1456	RSZ-1456
605 cc	15,0 cm	13,5 cm	6,6 cm	334-1506	RSZ-1506

Note: Individual implant dimensions may vary slightly in products of this type. Not all units will conform exactly to the dimensions noted above.
*Resterilizable up to 10 times.

INTERNATIONAL OPERATIONS



201 Mentor Drive
Santa Barbara, CA 93111 USA
Tel: +1 805 879 6000
Fax: +1 805 967 7108
www.mentorcorp.com

Manufacturer:
Mentor Medical Systems B.V.
Zernikedreef 2
2333 CL Leiden
The Netherlands

Mentor Medical Systems, Australia, Pty. Ltd. Tel: +61 2 9251 8111 Toll Free: +1 800 338 160 Fax: +61 2 9251 8333
Mentor Medical Systems, Canada, Inc. Tel: +1 905 725 7763 Toll Free: +1 800 668 6069 Fax: +1 905 725 7340
Mentor Deutschland GmbH Tel: +49 811 600 50 0 Fax: +49 811 600 50 11
Mentor Medical Systems Iberica, S.L. Tel: +34 91 562 2700 Fax: +34 91 728 3367
Mentor Medical Italia, S.r.l. Tel: +39 0 2 880 7761 Fax: +39 0 2 7209 3710
Mentor Medical Systems Ltd., U.K. Tel: +44 1635 511800 Fax: +44 1635 511810
Peroue Plastie, S.A.S., France Tel: +33 3 44 08 26 26 Fax: +33 3 44 08 26 01